



Collaborative Study of the Charm ROSA Safe Level Aflatoxin M1 Quantitative (SLAFMQ) Lateral Flow Test for Raw Bovine Milk

Robert S. Salter, David Douglas, Mark Tess, Robert J. Markovsky and Steven J. Saul, Charm Sciences, Inc., Lawrence, MA, USA 01843 bobs@charm.com

Introduction:

A 21 laboratory collaborative study tested raw bovine milk containing aflatoxin M1 using the Charm ROSA Safe Level Aflatoxin M1 Quantitative (SLAFMQ) lateral flow method. A ROSA reader quantitatively interpreted test strips with ppt readings. Readings lower than or equal to 400 ppt were negative, and readings greater than 400 ppt were interpreted positive. Initial positive samples were subsequently assayed two times. If both retests were greater than 400 ppt, the sample was called positive/actionable at 500 ppt (1,2).

Method:

Aflatoxin M1 standard (Sigma # 49319-U) used to prepare raw milk at 300, 350, 400, 450, 500 and 550 ppt.

Duplicate samples were blind coded, shipped on ice and tested within 1 week. Data were reported to a third party unaffiliated with sample preparer, Charm Sciences Inc

Laboratories performed the SLAFMQ method as follows in Figure 1

Figure1: Summary of SLAFMQ method

1. Place strip into ROSA incubator at 56°C.

2. Dilute raw milk into cold dilution buffer, mix and pipet 300 µL into strip. Seal Strip and close incubator lid which starts an 8 min timer.

3. After 8 minutes, place strip into ROSA reader for 5 sec analysis. The reader displays the determined part per trillion concentration. Values greater than 400 ppt are positive and need retesting in duplicate. Two more duplicate results greater than 400 ppt indicate an actionable sample.




Figure 2 Legend: Laboratories quantified results are plotted on x-axis versus prepared concentration on y-axis. Graph is divided into 4 quadrants by the blue vertical 400 ppt reader control point line and by the dotted horizontal 500 ppt action level line. Data points are true positive (TP) if they appear in red area, true negative (TN) if they appear in green area, false negative (FN) if they appear in blue area, and false violative-containing aflatoxin less than the action level (FP) if they appear in yellow area. The red curve plots the initial result percent positive (# positive/#tested x 100) on secondary y-axis versus concentration.

Discussion:

The Charm SL Aflatoxin M1 Quantitative (SLAFMQ) assay for raw milk in a 21 laboratory collaborative study detected US and Codex action levels with a 90% positive concentration with 95% confidence at 503 ppt and a 4.8% false negative rate. The initial test results prior to retesting have a 470 ppt 90/95 detection level with 2% false negative rate, which meet Interstate Milk Shipments single laboratory method parameters (7). False violative results were reduced to half from the initial assay result using a confirmation procedure that required samples with initial positive results to be retested twice, and for both retests to be positive.

The calculated repeatability (*r*) and reproducibility (*R*) for the SLAFMQ method were lower than published values for HPLC methods at comparable concentrations (8). These statistical parameters represent the 95% variation range expected from identical sample determinations within a laboratory (*r*) and between different laboratories (*R*). Lower *r* and *R* values and very low Horrat values, Horrat_r < 0.3 and Horrat_R < 0.5, indicate the SLAFMQ method has precision equivalent or better to HPLC for quantitation of aflatoxin M1 levels in milk. The SLAFMQ method had greater confidence at detecting actionable samples, at 500 and 550 ppt, than HPLC methods performed on the same samples.

Results:

The initial first sample results are plotted in Figure 2. The combined laboratory results of the SLAFMQ method, including positive retesting, and the initial first test result are presented in Table 1. The 90% positive level with 95% confidence is calculated from probit analysis (3). The statistical parameters of the initial test SLAFMQ quantitative reading are presented in Table 2 (4).

HPLC analyses of study samples were performed by five laboratories with one of those laboratories analyzing the samples twice (5, 6). A 38% false negative rate was found since 9 of 24 samples at 500 or 550 ppt aflatoxin M1 were reported to contain less than the 500 ppt action level. There were no false violative results in the final reported HPLC determinations.

Table 1: Dose Responses of SLAFMQ Method and Initial Test with 90% Positive Concentrations with 95% Confidence

Added Concentration (ppt)	# Samples Tested	SLAFMQ Method		Initial Test	
		# positive	% positive	# positive	% positive
0	41 ^A	0	0%	0	0%
300	42	0	0%	2	5%
350	42	9	21%	22	52%
400	42	6	14%	21	50%
450	42	39	93%	42	100%
500	42	39	93%	40	95%
550	41 ^A	40	98%	41	100%
90% Positive Concentration with 95% Confidence (one tail)		503 ppt		470 ppt	
Pearson Chi Square		25.7		20.2	

^AOutliers at 0 and 550 removed

Table 2: Charm SL Aflatoxin M1 Quantitative Test (SLAFMQ) statistical parameters showing mean, repeatability, reproducibility, and Horrat statistics

Added conc. (ppt)	Mean (#) of SLAFMQ Determinations	Intra-laboratory Repeatability Statistics				Inter-laboratory Reproducibility Statistics			
		StdDev _r	CV, % (RSV _r)	Repeatability <i>r</i> = 2.8 (StdDev _r)	Horrat _r Value (RSV _r /PRSV _r)	StdDev _R	CV _R % (RSV _R)	Reproducibility <i>R</i> = 2.8 (StdDev _R)	Horrat _R Value (RSV _R /PRSV _R)
0	6	6	119%	16	1.25	13	282%	36	2.73
300	291	31	11%	98	0.20	41	14%	115	0.28
350	388	61	16%	170	0.30	78	20%	218	0.39
400	394	40	10%	112	0.20	50	13%	139	0.25
450	505	48	9%	133	0.19	63	12%	176	0.25
500	495	62	13%	173	0.25	82	13%	173	0.25
550	595	35	6%	101	0.12	45	8%	125	0.15

References:

- Office of Regulatory Affairs, Section 527.4 CPG 7106.10, 11/29/05 http://www.fda.gov/ora/compliance_ref/cpg/cpgfdg/cpg527-400.html
- Codex Alimentarius Commission, Standard 232-2001 Aflatoxin Levels in Milk, http://www.codexalimentarius.net/download/standards/399/CXS_232e.pdf
- XL-Stat, Addinsoft, 224 Centre Street, 3rd Floor, New York, NY 10013, USA Tel: 866 740 STAT (toll free), Fax: (646) 349 1715 Web: www.xlstat.com
- International Standard Organization, (1994) International Standard Organization, Brussels Belgium, 5725-2
- Official Methods of Analysis of the AOAC INTERNATIONAL (2005) 18th Ed., AOAC INTERNATIONAL, Gaithersburg, MD, 2000.08
- Official Methods of Analysis of the AOAC INTERNATIONAL (2005) 18th Ed., AOAC INTERNATIONAL, Gaithersburg, MD, 986.16
- Center for Veterinary Medicine Data Requirements for Milk Screening Tests Labeled for Testing Milk Tankers at the Bulk/Tank Tanker Truck for Drug Residues Tests with Instrument Readers/Printers Only- January 31, 2002
- Tuinstra, L.G.M.T., Roos, A.H., and van Trijp, J.M.P. (1993) Journal of the AOAC. 76, 1248-125

Acknowledgment: The authors thank the collaborators and Mr. Gene Wright of DSHS- Austin Texas for his assistance in collecting, scoring and decoding data.

Collaborators: C. Hall and M. Kingsley of Texas DSHS-Austin; G. Dixon and N. Turnage of Tarrant County Public Health; K. McClaran and K. Swacina of City of Lubbock Public Health; M. Williams and S. Wolff of San Antonio Metro Health Dept.; C. Lucero and C. Fierro of El Paso City Env. Health District; K. Keys of Northeast Texas Public Health District; A. Brandt of Pioneer Labs-Artesia; P. Gutierrez and A. Reeb of NM Dept. of Agr.; K. Jones and R. McDaniel of NC Dept. of Agr. Food and Drug Protection; T. Tieso of NE Dept. of Agr.; K. Challquist of KS Dept. of Agr.; P. Stetzel of Iowa Dept. of Agr.; A. Gilleland of GA Dept. of Agr.; D. Klein of AZ Dept. of Agr.; D. Manor of AR Dept. of Agr.; F. Barcellos of OK Dept. of Agr.; B. Malone of Trilogy Laboratory; B. Shaffer of DFA- Kansas City; T. Hickerson and T. Moberly of DFA-Graprevine; R. Lam and J. Miller of FL Dept. or Agr.; J. Nandrea, of U.S.F.D.A. Kansas City District Laboratory