Q-194





Visit Us at Booth #1106

Purpose

The purpose of this study is to simplify coliphage detection to meet drinking water standards within one working day. Described are a presence-absence (PA) method for 100 ml samples, MPN quantitation techniques for 1-100 ml samples, and a continuous detection technique over 24 hours, testing liter volumes.

Background

Recently, coliphage – viruses to the coliform group of bacteria – were listed as an equivalent fecal marker to E. coli in the Ground Water Rule (GWR). Half of waterborne illnesses are viral in origin. Public health officials agree that bacterial indicators and viral indicators provide valuable information about water sources and processes. Comprehensive risk assessment of water should include tests for both bacterial and viral indicators. Conventional coliphage detection methods have been difficult, time consuming and expensive. Bacterial indicator tests have been used predominantly because of ease of use and low cost. Simple, fast and low-cost coliphage methods are needed to compliment bacterial testing.

Same Day and Continuous Detection of Coliphage in Water Robert S. Salter, Gregory W. Durbin, Shishan Wang, Rami Abraham, and Stanley Charm Charm Sciences, Inc., 659 Andover St., Lawrence, MA USA Phone 800-343-2170 Fax 978-687-9216 email: bobs@charm.com

Materials and Methods

PA Method: Modifications to Method 1601, called Fast Phage™, are summarized in Table 1. Wastewater-spiked water samples at 0.5 pfu/100 ml were split into 20 replicates and tested by Method 1601 and Fast Phage. After enrichment step (5hrs) a transfer to MUG-gal media predicts by fluorescence (366 nm) the next-day plaque results in less than 3 hours (total time from assay start <8 hours). A variety of secondary wastewaters pre-and post-chlorine disinfection were used to spike ground water.

MPN Methods – TEMPO[®] Method: 1 to 3.7 ml water sample plus 300 ul of 5 ml saline rehydrated CN-13 *E. coli* tablet to final 4 ml volume are added to TEMPO[®] card. **100 ml Method:** 1-95 ml sample to final 95 ml volume and add STEP-2 pouch. Sample is pre-heated in 38°C water bath while E. coli tablet is rehydrated for 10 minutes in 5 ml saline. Contents are combined and sample divided. In this experiment, Quanti-Tray[®] MPN are performed. Comparative Methods 1601, 1602 or Double Layer Agar are also performed on same samples. MPN devices are incubated 5.5 - 6 hours. MPN are determined from fluorescent (366 nm) sub-divisions and using MPN table.

Continuous Method: A continuous monitor of coliphage was constructed, Figure 1, by culturing E. coli host in a continuously fed water sample/media mix in a disposable container with external monitoring probes. Transmittance 520 nm 20 cm path length T=0.25 is measured by a probe/reflector. Software monitors transmittance and activates peristaltic pump that in-feeds mix and out-feeds E. coli culture to maintain constant T Fluorescence 366 nm (450 nm excite) is monitored by Turner Designs C-7 probe. 1.0 mV activates warning of coliphage host infection. A mix pump combines sample with 10x media concentrate and a 90°C heated coil is used for culture disposal.

Table 1: Summary of Fast Phage[™] Modifications Compared to Method 1601

| METHOD 1601 | FAST PHAGE MODIFICATION |
|--|---|
| Pre-culture <i>E. coli</i> CN-13 overnight. Transfer to TSB and bring to OD 520 _{nm} = 0.1 - 0.5. Ice. | Pre-dispensed dried media components in dissolvable film. Formulation proprietary for rapid |
| Prepare 10x TSB, MgCl ₂ , Naladixic Acid solutions in advance and autoclave. | amplification. Film pouch also contains <i>E. coli</i> tablet (CN-13 for somatic and F_{amp} for F+) inoculum that dissolves. |
| Add prepared solutions and <i>E. coli</i> to water sample in precise sequence of additions and temperatures. | Add media/ <i>E</i> . <i>coli</i> pouch to water sample in test vessel. Gently swirl to dissolve contents ~ 10 minutes. |
| Incubate overnight at 36° ± 1°C air incubator. | Incubate 30 minutes in water bath at $38^{\circ} \pm 1^{\circ}$ C. Incubate 4 hours and 30 minutes in air incubator at $39^{\circ} \pm 1^{\circ}$ C. |
| Use pre-cultured log phase <i>E. coli</i> to make plaque plates. Prepare TS agar, autoclave and add pre-prepared and sterile filtered antibiotic solution after temperate. | Use <i>E. coli</i> tablet as TSB inoculum. Use supplied freeze- dried antibiotic (naladixic acid for somatic or ampicillin-streptomycin for F+) to make solution. Prepare TS agar, autoclave and add antibiotic solution after temperate. |
| Transfer 10 ul of enrichment culture (24 hours) to plaque plate. Incubate overnight at 36° ± 1°C air incubator. | Transfer 10 ul of enrichment culture (5 hours) to plaque plate. Incubate overnight at 39° ± 1°C air incubator. |
| NO RAPID PREDICTION OF PLAQUE RESULT. | Transfer enrichment culture (1 ml for F+ or 10 ml somatic) to 100 ml indicating medium STEP-2. Incubate up to 3 hours for fluorescence positive to predict positive plaque result. |

Figure 1: Continuous Testing System




Figure 1 Legend: 60 - 150 ml bacterial culture in a disposable flask in probe holder on a heated base plate at 38°C. Media/water-sample mix is fed in and culture removed at the same rate using disposable tubing and external transmittance probe-reflector with computer software (not shown). External fluorometer probe detects coliphage infection of culture.

Results & Discussion

Table 2 compares positives of 20
 replicates detected by Method 1601 to those detected by Fast Phage STEP-1 plaque positive and fluorescence.

Table 3 relates the accuracy of the
 fluorescent indicator at predicting the STEP-1 plaque results. The predictive correlation is 90%. Disinfected waters had 7 false negatives for fluorescence, while untreated waters had 5 fluorescence positives without STEP-1 plaque confirmation, but with detected coliphage from STEP-2 plaque.

Table 4 reports somatic phage TEMPO
 MPN determinations of non-disinfected secondary waste waters compared to Method 1602 of plaque forming units (pfu).
Table 5 reports somatic phage by 100 ml
 MPN technique. Results are in agreement within 0.5 log. CVs of replicate MPN are about 30 - 50%.

Figure 2 demonstrates a continuous overnight sampling of over 10 L water with a next morning challenge of 20 pfu F+ coliphage in 1 L using F_{am} E. coli host. Positive fluorescence is 5.5 hours post-phage challenge. Plaque were confirmed from fluorescent positive sample. Similar results using CN-13 host are observed with a daily flow rate of about 4-6 liters/day. No fluorescence is observed in unchallenged water samples with up to 20 L sample passing through the culture flask over the 24 - 48 hour sampling period. Biofouling at 24 - 48 hours requires disposal and exchange of flask/tubing.

| Table2: Comparability Study | | | | | | | |
|-----------------------------|---------------|-------------|-----------------|----------|--------|--------------|---|
| | Waste | Chlorine | Dilution | EPA 1601 | STEP-1 | Fluorescence | |
| | Location | Disinfected | | Plaque | Plaque | Positive | S |
| | M 10-15-08 | no | 2.25 ml / 4.5 L | . 15 | 15 | 16 | |
| | L 09-09-08 | no | 4.5 ml / 4.5 L | . 12 | 10 | 11 | |
| | SA 12-18-08 | no | 0.35 ml / 4.5 L | . 10 | 6 | 7 | |
| | L 06-26-08 | yes | 10 ml / 4.5 L | . 15 | 13 | 7 | |
| | O 12-30-08 | no | 1.9 ml / 4.5 L | . 10 | 9 | 11 | |
| | O 01-02-09 | yes | 280 ml / 4220 L | . 5 | 4 | 3 | |
| | Sum of Positi | Ves | | 67 | 57 | 55 | |
| | % Positive | | | 55.8% | 47.5% | 45.8% | |

| | Table J. Fludiescence Fledicilon Di STEF-I Flaque Results | | | | | | | |
|---|---|---------------|----------|---------|---|---|--|--|
| | Fluorescence | STEP-1 Plaque | | Total E | | Ρ | | |
| ľ | | Positive | Negative | Total F | Sensitivity | 8 | | |
| | Positive | 50 | 5 | 55 | Specificity | Q | | |
| | Negative | 7 | 58 | 65 | False Positive Error | | | |
| | Total STEP-1 | 57 | 63 | 120 | False Negative Error Overall Agreement | ç | | |


Table 4: MPN Determination of Somatic Phage in ml of Secondary Wastewaters
 ompared to 1602 Plaque Determination

| Secondary Wastewater | Method 1602 | | TEMPO | | | |
|----------------------|---------------|-------|---------------|-----|--|--|
| Sample 1 ml/test | Average (N=4) | Stdev | Average (N=6) | Sto | | |
| M 1015 | 52 | 13 | 42 | 18 | | |
| L 0909 | 12 | 4 | 10 | 4 | | |
| 10623 | 34 | 8 | 39 | 1 | | |
| l 0623 1:5 | 6.5 | 2 | 6 | 2 | | |
| O 1230 | 16 | 3 | 14 | 8 | | |
| V 1223 | 56 | 8 | 65 | 4 | | |
| L 012209 | 526 | 56 | 197 | 6 | | |
| x174/ml PC1118 | 77 | 19 | 49 | 1 | | |
| x174/bottle PC 1111 | 173 | 11 | 101 | 28 | | |

Table 5: MPN Determination of Somatic Phage in 1 ml of Secondary Wastewaters mpared to Double Layer Agar and Method 1601 MPN

| Waste Water | Quanti-Tray MPN/mI Average N=2 | Double Layer pfu/ml | Enumera MPN 1 pfu/n | |
|-----------------|--------------------------------------|---------------------------|---------------------------|--|
| L 01-23-08 | 163.6 | 107 | 66.7 | |
| L 01-23-08 1:10 | 14.5 | 10.7 | 6.7 | |
| W 02-01-08 | 4.1 | 12.5 | 4.5 | |
| M 01-15-08 | 26.7 | 40 | 10 | |

ransmittance Signal of F_{am} Continuous



Sample Number (Minutes)




Figure 2a Legend: 60 ml culture of *E. coli* F_{amp} was grown to ansmittance 0.25 at which time media xed with water sample was fed into the overnight. The next morning (sample 1280) 20 pfu F+ coliphage in 1 liter was introduced was reconnected


Figure 2b Legend: Fluorescent reading of the *E. coli* during their growth and continuous state. hage infection starts to be detected a sample (minute) 1687 and by sample warning signal trips that infection ha

Figure 2c Legend: Flow rate is 0 ml/min until transmittanc meets target of 0.25 when pump is vivated Rate is determined by softwar Igorithm that has a high-low flow inflection at the target transmission. As transmission gets less than target, flow jumps. When ansmission is higher than target, flow rate





Summary







Fecal assessment by detecting less than 1 pfu coliphage in 100 ml water meets ground water detection guidelines and can now be done in less than 8 hours.

Fast Phage 100 ml P/A modifications simplify Method 1601 to a 16-hour procedure with an 8-hour 90% correlative fluorescent prediction.

Somatic MPN techniques offer 6-hour fecal indicator assessment of fresh, untreated source water, providing same-day pollution-risk assessment.

Continuous sample monitoring for coliphage is achievable with 4-6 hour detection. This offers liter volume testing of continuously sampled water as well as automation of fecal risk assessment.

References

- National Primary Drinking Water Regulations: Ground Water Rule [Federal Register: November 21, 2006 (Volume 71, Number 224)]
- 2. Method 1601: Male-specific (F+) and Somatic Coliphage in Water by Two-step Enrichment Procedure, EPA Number 821R1030 April, 2001
- . Method 1602 Male-Specific (F+) and Somatic Coliphage in Water by Single Agar Layer (SAL) Procedure, EPA Number 821R1029 April, 2001